Medix Biochemica

Anti-hCG beta 5006 SP-5

Product overview									
Catalog number	100005								
Specificity	Antibody recognizes human chorionic gonadotropin and its free beta subunit								
Description	Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components.								
Product buffer solution	0.9 % NaCl, 0.095 % NaN $_3$ as a preservative								
Shelf life and storage	24 months from manufacturing at 2–8 °C								
Subclass	IgG ₁								
Analyte description	Human chorionic gonadotropin (hCG) is a glycoprotein hormone product in pregnancy by the developing embryo soon after conception and later the syncytiotrophoblast (part of the placenta). Its role is to prevent the disintegration of the corpus luteum of the ovary and thereby maintain progesterone production that is critical for a pregnancy in humans. Earl pregnancy testing, in general, is based on the detection of hCG. hCG is produced also by some tumors, but it is not known whether this produce is a contributing cause or an effect of tumorigenesis.								
Parameters tested on each lot									
Product appearance	Liquid, may turn slightly opaque during storage								
Product concentration	5.0 mg/ml (+/- 10%)								
Immunoreactivity	80–120% compared to the reference sample in an FIA test								
IEF Profile	6.2–7.4								
Purity	≥ 95 %								
Kinetic parameters									
Association rate constant	hCG: 5.7 x 10 ⁶ 1/Ms and hCGβ: 2.3 x 10 ⁶ 1/Ms								
Dissociation rate constant	hCG: 2.9 x 10 ⁻⁴ 1/s and hCGβ: 4.8 x 10 ⁻⁴ 1/s								
Affinity constant	hCG: K _A = 2.0 x 10 ¹⁰ 1/M; K _D = 5.0 x 10 ⁻¹¹ M (= 0.05 nM) hCGβ: K _A = 4.8 x 10 ⁹ 1/M; K _D = 2.1 x 10 ⁻¹⁰ M (= 0.21 nM)								
Determination method	SPR analysis (ProteOn XPR36)								
Determination antigen	hCG, Scripps (Cat C0714, Lot 2430801); hCGβ, Scripps (Cat C0914, Lot 2310001)								



Oy Medix Biochemica Ab Klovinpellontie 3, FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

Medix Biochemica

Cross-reactivities

LH (K_D = 24 nM) Does not recognize hCG α , FSH, or TSH

Epitope

Beta-2 as described in Berger et al. (2013). The antibody recognizes both intact hCG and free β subunit.

DETECTION

Pair recommendations

						hCG	beta				alp	alpha subunit	
			5004	5006	5008	5009	5011	5012 free β	5014	5016	5501	5503	6601
		5004	-	-	-	+	+	-	+	-	+	+	+
		5006	-	-	-	-	-	-	+	-	+	+	+
	CAPTURE	5008	-	-	-	+	-	-	+	-	+	+	+
		5009	+	+	+	-	-	-	+	+	-	-	+
		5011	+	+	+	-	-	-	+	+	-	-	+
		5012 free β	+	+	+	-	-	-	+	+	-	-	-
		5014	+	+	+	+	+	-	-	+	+	+	+
		5016	-	-	-	-	-	-	+	-	+	+	+
	Please note that pair recommendations are based on results obtained by our laboratory. Equally good results may be obtained using other pairs and therefore these recommendations are only indicative.												
Platforms tested	FIΑ,	CLIA											
Antigens tested	Native hCG antigens, Lee Biosolutions 189-10 and 189-11 Native β-hCG antigen, Lee Biosolutions 325-11												
Product stability	TEMPERATURE, TIMERESULTNot Determined (N/D)-												
Miscellaneous	-												
References	Berger, P., Paus, E., Hemken, P.M., Sturgeon, C., Stewart, W.W., Skinner, J.P., Harwick, L.C., Saldana, S.C., Ramsay, C.S., Rupprecht, K.R., Olsen, K.H., Bidart, J.M. and Stenman, U.H. (2013) Candidate epitopes for measurement of hCG and related molecules: the second ISOBM TD-7 workshop. Tumor Biol., 34: 4033-4057.												



Oy Medix Biochemica Ab Klovinpellontie 3, FI-02180 Espoo, Finland medix@medixbiochemica.com www.medixbiochemica.com VAT reg. no. FI14631532

Legal disclaimer

Medix Biochemica products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.