

Anti-Influenza B 9906 SPTN-5

Product overview

| | |
|--------------------------------|---|
| Catalog number | 100437 |
| Specificity | Antibody recognizes <i>Influenza B</i> nucleoprotein (np) |
| Description | Monoclonal mouse antibody, cultured <i>in vitro</i> under conditions free from animal-derived components. |
| Product buffer solution | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| Shelf life and storage | 12 months from manufacturing at 2–8 °C |
| Subclass | IgG _{2b} |
| Analyte description | Influenza B virus is the only species in a genus in the virus family Orthomyxoviridae. Influenza B viruses are only known to infect humans and seals. This limited host range is apparently responsible for the lack of Influenza B caused influenza pandemics in contrast with those caused by the morphologically similar Influenza A. In addition, Influenza B mutates at a rate 2-3 times lower than type A. However, influenza B mutates enough that lasting immunity is not possible. |

Parameters tested on each lot

| | |
|------------------------------|--|
| Product appearance | Liquid, may turn slightly opaque during storage |
| Product concentration | 5.0 mg/ml (+/- 10 %) |
| Immunoreactivity | 80–120 % compared to the reference sample in an FIA test |
| IEF Profile | 6.5–7.9 |
| Purity | ≥ 95 % |

Kinetic parameters

| | |
|-----------------------------------|---|
| Association rate constant | 7.6×10^4 1/Ms |
| Dissociation rate constant | 8.0×10^{-4} 1/s |
| Affinity constant | $K_A = 9.5 \times 10^7$ 1/M; $K_D = 1.1 \times 10^{-8}$ M (=11.0 nM) |
| Determination method | SPR analysis (ProteOn XPR36) |
| Determination antigen | Recombinant Influenza B virus nucleoprotein from strain B/Singapore/222/79 (in-house antigen, UniProtKB P04666) |



Legal disclaimer

| Cross-reactivities | Recognizes Influenza B/Yamagata clade 3, Influenza B/Victoria clade 1A, clade 1A ($\Delta 2$), and clade 1A ($\Delta 3$). Does not recognize Influenza A (H3N2) group 3C.2a1b and group 3C.3a, Influenza A (H1N1)pdm09 group 6B.1A1 and group 6B.1A5. | | | | | | | | | | | | |
|-----------------------------|---|-------------------|--------|-----------------|----|-----------------|----|----------------|----|-----------------|----|----------------|----|
| Epitope | Not determined (N/D) | | | | | | | | | | | | |
| Pair recommendations | As nucleoprotein is abundant in the virus, all Anti-Influenza B antibodies (9901, 9906, and 9909) pair with themselves and each other. | | | | | | | | | | | | |
| Platforms tested | FIA | | | | | | | | | | | | |
| Antigens tested | Recombinant Influenza B NP antigen, Medix Biochemica 610051 | | | | | | | | | | | | |
| Product stability | <table><thead><tr><th>TEMPERATURE, TIME</th><th>RESULT</th></tr></thead><tbody><tr><td>-70 °C, 21 days</td><td>OK</td></tr><tr><td>-20 °C, 21 days</td><td>OK</td></tr><tr><td>+4 °C, 21 days</td><td>OK</td></tr><tr><td>+35 °C, 21 days</td><td>OK</td></tr><tr><td>+45 °C, 3 days</td><td>OK</td></tr></tbody></table> <p>Stability testing is performed in the product buffer to see whether different temperatures affect the antigen binding, charge or composition of the antibody. Please note that the shelf life given on the first page is based on real time stability testing at 2–8 °C in the product buffer.</p> | TEMPERATURE, TIME | RESULT | -70 °C, 21 days | OK | -20 °C, 21 days | OK | +4 °C, 21 days | OK | +35 °C, 21 days | OK | +45 °C, 3 days | OK |
| TEMPERATURE, TIME | RESULT | | | | | | | | | | | | |
| -70 °C, 21 days | OK | | | | | | | | | | | | |
| -20 °C, 21 days | OK | | | | | | | | | | | | |
| +4 °C, 21 days | OK | | | | | | | | | | | | |
| +35 °C, 21 days | OK | | | | | | | | | | | | |
| +45 °C, 3 days | OK | | | | | | | | | | | | |
| Miscellaneous | - | | | | | | | | | | | | |
| References | - | | | | | | | | | | | | |



Legal disclaimer